

We have made further advancements for you

and are pleased to present to you the:

## IMPROVED SALIVA DRUG TEST RapidSTAT®:

New incubation device:

**Easier closure of the lid**

Larger openings of the incubation device:

**Easier and improved pipeting of the test**

Darker negative lines:


**Negative lines already visible from the 3. minute on  
And easy to read in the 8. minute.**

Positive results stable for 20 minutes:

**Test can be evaluated from the 6.-8. minute -  
Stable results up to 20 minutes !**

**AVAILABLE  
NOW !**

The new RapidSTAT®  
is available since May 2011



**MAVAND**  
SOLUTIONS GMBH

Ulrichstr. 21  
72116 Mössingen  
Germany

fon +49(0)7473-95 80-28  
fax +49(0)7473-95 80-29  
email info@mavand.com  
internet www.mavand.com