



RAPID Saliva SET™

The **RAPID Saliva SET™** in a **Clip-System** is a new system designed for the collection, storage and transport of saliva samples in connection with the detection of illegal drugs by standardised laboratory analysis (immuno assay, GC-MS, LC-MS).



Technology

The patent pending aroma label ensures that in case of less saliva of the donor the salivation will be activated by the scent.

The new developed patent pending sample collector is furnished with a unique fine micro fiber which allows to collect quickly the oral fluid. It is not necessary to collect a lot of saliva, only 0,15 up to 0,3 ml are sufficient to find all possible concentrations.

The buffer solution assures that the fat binding drug will be washed out of the micro fiber. No extensive procedure is necessary for the extraction. The buffer saliva mixture can be used directly.

User friendly

The **RAPID Saliva SET™** in a **Clip-System** can be used easily, quickly and hygienically in 3 steps. Collect the sample for 30 seconds, wash out the swab for 10 seconds, put the stick in its tube and close it. The collection can be observed easily. There is no invasion of privacy of the donor.

Application field

The **RAPID Saliva SET™** in a **Clip-System** can be applied for drug monitoring during prevention campaigns, rehabilitation measures, access controls and security examinations as well as in case of suspicion of substance abuse on work places.

Service

The sample can be sent to the analytical service of MAVAND Solutions GmbH or to another commissioned laboratory using its primary packing.

