



Rapid STAT®

Rapid STAT® is a **revolutionary point of collection test** for the detection of drugs of abuse in oral fluid, on surfaces or for pure substance testings.

Rapid STAT® combines the **convenience** of oral fluid collection, surface wipe testing or pure substance measurements with the **sensitivity, accuracy** and **precision** of a laboratory based test providing results in **about 6-13 minutes**.

Patented method



AMP

Amphetamines 25 ng/ml

BZO

Benzodiazepines 25 ng/ml

COC

Cocaine 12 ng/ml

MTD

Methadone 15 ng/ml

MET•MDMA

Methamphetamines 25ng/ml•
Ecstasy 50 ng/ml

OPI

Opiates 25 ng/ml

THC

Marijuana 15 ng/ml

FLEXIBILITY

- 7 drugs available
- 2 to 7 drugs per test - varied configurations
- Can be applied as saliva-, wipe- or pure substance test

CONVENIENCE

- Easy to use due to the One-Hand-Clip-System
- All collections can be observed
- No restroom needed
- Specimen collection takes just 30 seconds
- Increased salivation due to innovative aroma field
- No shy bladder issues
- No specimen dilution or tampering issues
- **Results in about 6- 13 minutes**
- Rapid results allow for correct decisions

PERFORMANCE

- Each test strip has built-in procedural control
- Colloidal gold technology = superior sensitivity, accuracy & precision
- Tested for cross-reactivity against most common drugs and other substances
- Detects parent THC

SAMHSA states: "The only way to detect marijuana use is through the presence of the parent drug (THC) in the oral fluid....."

TRAINING/SUPPORT

- Onsite training
- Video Training on DVD
- Telephonic technical support