



NEW

Rapid WIPE®

Rapid WIPE® convinces due to the easy handling, sensitivity and accuracy. The results are available in 3-5 minutes !



AMP

Amphetamines

BZO

Benzodiazepines

COC

Cocaine

MET

Methamphetamines

MDMA

Ecstasy

OPI

Opiates

THC


Marijuana

FLEXIBILITY

- Detection of up to 6 drugs
- 1 to 6 drug parameter per test - several configurations
- Can be used as surface wipe test or substance test

USER FRIENDLY

- Easy to handle - ready in just 4 steps
- Sample collection in just 10 seconds
- Flexible, long collection stick
- **Results in just 3-5 minutes**
- Easy and save interpretation of results
- Results are stable for 20 minutes and longer
- Short instruction on each single packing
- Small handy package (LxB 23 cm x 8 cm)



CHARACTERISTICS

- Each test strip has an integrated function control
- Tested on cross reactivity of the most common drugs and other substances
- Detects THC direct (not metabolite)

TRAINING / SERVICE

- Training onsite
- Technical telephone support