

NEWSLETTER

Ausgabe: 0032

LIGHTBOX for MAVAND test kits!

MAVAND designs special tool to facilitate the interpretation of test results for emergency task forces!


The drugs of abuse test kits Rapid STAT[®], Rapid WIPE[®] and Rapid TOX[®] are used by police task forces, military, customs and security forces around the world. Very often, poor lighting conditions avoid the efficient and reliable interpretation of the test results. For this reason, MAVAND designed the LIGHTBOX which is equipped with LEDs and a special light filter. The green light enables a better interpretation of the negative lines of the Rapid STAT[®], Rapid WIPE[®] and the Rapid TOX[®] in the dark or if lighting conditions are bad.



LIGHTBOX designed for mobile use:

- (BxHxD): 70x125x25mm
- Weight: 0.14kg
- Battery-powered: AA Batteries, 3.0V / 0.02A
- Max. 0.06 Watt - Interpretation of at least 5'000 test results x approx. 20 sec. use
- Storage and operating temperature -20°C to 50°C (please consider the Rapid STAT[®], Rapid WIPE[®] and Rapid TOX[®] operating temperatures)

Applicable for Rapid STAT[®], Rapid WIPE[®] and Rapid TOX[®]



Please do not hesitate to contact us for free test trials!



MAVAND
SOLUTIONS GMBH

Ulrichstr. 21
72116 Mössingen
Germany

fon +49(0)7473-95 80-28
fax +49(0)7473-95 80-29
email info@mavand.de
internet www.mavand.de